



# Tower Crane

USB2.0 version

## *User's Manual*



[www.inteco.com.pl](http://www.inteco.com.pl)

ver. 9.6

---

## **COPYRIGHT NOTICE**

---

**© Inteco Limited**

All rights reserved. No part of this publication may be reproduced, stored in a retrieval system, or transmitted, in any form or by any means, electronic, mechanical, photocopying, recording or otherwise, without the prior permission of Inteco Ltd.

---

## **ACKNOWLEDGEMENTS**

---

Inteco Ltd acknowledges all trademarks.

IBM, IBM - PC are registered trademarks of International Business Machines.

MICROSOFT, WINDOWS 2000/XP/7 are registered trademarks of Microsoft Corporation.

MATLAB, Simulink and RTW are registered trademarks of Mathworks Inc.

## Contents

<b>1. INTRODUCTION AND GENERAL DESCRIPTION.....</b>	<b>5</b>
1.1. Product overview.....	6
1.2. Requirements.....	8
<b>2. STARTING, TESTING AND STOPPING PROCEDURES.....</b>	<b>9</b>
2.1. Starting procedure.....	9
2.2. Testing and troubleshooting.....	9
2.3. Stopping procedure.....	14
<b>3. MAIN CONTROL WINDOW.....</b>	<b>15</b>
3.1. Tools.....	15
3.2. Drivers.....	22
3.3. Simulation Model.....	22
3.4. Demo Controllers.....	23
<b>4. YOUR FIRST REAL-TIME CONTROL EXPERIMENT.....</b>	<b>26</b>
4.1. Real-time experiment.....	31
4.2. Data processing.....	32
<b>5. SIMULATION .....</b>	<b>34</b>
5.1. Simulation control experiments.....	34
5.2. PID simulation control of load position.....	38
<b>6. PROTOTYPING YOUR OWN CONTROLLER IN THE REAL-TIME ENVIRONMENT.....</b>	<b>41</b>
6.1. Creating a model .....	43
6.2. Code generation and the build process.....	45
<b>7. MATHEMATICAL MODEL OF THE TOWER CRANE.....</b>	<b>49</b>
7.1. Equations.....	49
7.2. Comparison between the mathematical model and laboratory system.....	50
<b>8. DESCRIPTION OF THE TOWER CRANE CLASS PROPERTIES.....</b>	<b>53</b>
8.1. BaseAddress.....	54
8.2. BitstreamVersion.....	54
8.3. Encoder.....	55
8.4. PWM.....	55
8.5. PWMPrescaler.....	55
8.6. ResetEncoder.....	56
8.7. RailLimit.....	56
8.8. RailLimitFlag.....	56
8.9. RailLimitSwitch.....	57
8.10. ResetSwitchFlag.....	57
8.11. Therm.....	57
8.12. ThermFlag.....	57
8.13. Time.....	57
8.14. Quick reference table.....	58
<b>9. HOW TO FILL IN THE COMPILATION SETTINGS PAGE.....</b>	<b>59</b>



# Tower Crane

**The industrial tower crane model controlled from PC.  
The control goal: to track a trajectory and not to swing the load.**

## 1. Introduction and general description

Tower Crane is a nonlinear electromechanical system having a complex dynamic behaviour and creating challenging control problems. It is controlled from a PC. Therefore it is delivered with hardware and software which can be easily mounted and installed in a laboratory. You obtain the mechanical unit together with the power supply and interface to the PC and the dedicated digital board configured in the Xilinx® technology. The software operates under MS Windows® using MATLAB® and RTW toolbox package.

Besides the hardware and the related software you obtain *User's Manual*.

The manual:

- shows step-by-step how to design and generate your own real-time controller in the MATLAB®&Simulink® environment,
- contains the library of ready-to use real-time controllers,
- includes the set of preprogrammed experiments.




Fig. 1.1 The Tower Crane setup

The **Tower Crane setup** (Fig. 1.1) consists of a payload hanging on a pendulum-like lift-line wound by a motor mounted on a cart.

The payload is lifted and lowered in the  $z$  direction. Both the arm and the cart are capable of horizontal motion: the cart in the radial  $x$  direction along the arm and the arm in the rotary direction. The angular position of the arm is expressed by the  $\theta$  angle. The payload attached to the end of the lift-line can move freely in three dimensions. The Tower Crane is driven by three DC motors.

There are five measuring encoders measuring five state variables: the cart co-ordinates on the polar coordinates plane, the lift-line length, and two deviation angles of the payload. The encoders measure movements with a high resolution equal up to 4096 pulses per rotation (ppr). These encoders together with the specialised mechanical solution create a unique measurement unit. The deviation of the load is measured with a high accuracy equal to 0.0015 rad.

The power interface amplifies the control signals which are transmitted from the PC to the DC motors.

The PC equipped with the RT-DAC/USB multipurpose digital I/O board communicates with the power interface. The whole logic necessary to activate and read the encoder signals and to generate the appropriate sequence of pulses of PWM to control the DC motors is configured in the Xilinx<sup>®</sup> chip of the RT-DAC/USB board. All functions of the board are accessed from the Tower Crane toolbox which operates directly in the MATLAB<sup>®</sup>&Simulink<sup>®</sup> environment.

### KEY FEATURES of Tower Crane:

- Three-dimensional laboratory model of industrial crane.
- A highly nonlinear MIMO system.
- It can be easily installed.
- There are high-resolution sensors – unique 2D angle measuring unit.
- The set-up is fully integrated with MATLAB<sup>®</sup>&Simulink<sup>®</sup> and operates in real-time in MS Windows<sup>®</sup>.
- Real-time control algorithms can be rapidly prototyped. No C code programming is required.
- The software includes complete dynamic models.
- *User's Manual* contains the library of basic controllers and a number of pre-programmed experiments which familiarise the user with the system in a fast way.
- It is ideal for illustrating complex nonlinear control algorithms.

### 1.1. Product overview

The tower crane is delivered in partially mounted form.

#### SETUP COMPONENTS

##### hardware

- mechanical unit
- interface and Power Interface Unit
- RT-DAC/USB I/O board (the PWM control logic is stored in the XILINX chip)



Fig. 1.2 Cart and 2D angle measuring unit




Fig. 1.3 X axis drive

#### **software**

- TowerCrane Toolbox operating in MATLAB®&Simulink® environment

#### **e-manuals**

- *Installation Manual*
- *User's Manual*

## 1.2. Requirements

The following minimum configuration is required:

### Hardware:

- The Tower Crane System including the mechanical unit and Power Interface unit.
- Computer system based on INTEL or AMD processor.
- Specialised RT-DAC/USB2 I/O board.

### Software:

- Microsoft Windows XP/W7x86 and MATLAB 32 bit with Simulink and RTW toolboxes (not included),
- MS Visual C++ compiler,
- or the MS Visual Express compiler (free of charge),
- or the Open Watcom 1.9 compiler [www.openwatcom.org] free of charge

**or**

- Microsoft Windows XP/W7x64t and MATLAB 64 bit with Simulink and RTW toolboxes (not included),
- Microsoft Software Development Kit (SDK) 7.1 (free of charge compiler)



**Details of the required software are available at:**

**[http://www.inteco.com.pl/support/Software\\_requirements.pdf](http://www.inteco.com.pl/support/Software_requirements.pdf)**



**Neither MATLAB nor the compiler can be installed in the Program Files directory (name of the directory cannot include space).**

Manuals:

*Installation Manual*

*User's Manual*



**The experiments and corresponding to them measurements have been conducted by the use of the standard INTECO systems. Every new system manufactured and developed by INTECO can be slightly different to those standard devices. It explains why a user can obtain results that are not identical to these given in the manual.**



## 2. Starting, testing and stopping procedures

### 2.1. Starting procedure

Assemble the external USB board, Control Interface and Tower Crane together. Invoke MATLAB and type:

*tcr*

The *USB\_TCrane\_Main* window opens (see Fig.2.4). The pushbuttons indicate an action that executes callback routines when the user selects a menu item.




Fig.2.4 The *USB\_TCrane\_Main* window

The window contains testing tools, device driver, simulation model and demo applications. You can see a number of pushbuttons ready to use.



**In the case if the RT-DAC/USB I/O board is used only one software module can communicate with USB board. Remember do not open more than one Simulink model at time.**

### 2.2. Testing and troubleshooting

This section explains how to perform the tests. These tests enable to check the correctness of the mechanical assembling and wiring. The tests have to be performed obligatorily after assembling the system. They are also necessary in a case of an incorrect operation of the system. The tests can be helpful to find an eventual reason of the system fail. The tests have been designed to validate the existence and sequence of measurements and controls. They do not relate to accuracy of the signals.

The operational space of the tower crane is illustrated in Fig. 2.5.




Fig. 2.5 Operational space of Tower Crane

First, you have to be aware that all signals are transferred in a proper way. Eleven checking steps are applied.

- Double click the *Basic Tests* button. The following window appears (Fig. 2.6):



Fig. 2.6 Tower Crane Basic Tests window

The first step is to check the proper operation of the limit switches. There are three switches applied to stop the system motions and to secure the system against destruction if the cart or the arm approach the limits. The z axis motion switch is a typical mechanical one, two other switches are contactless activated if a metallic object comes closer to the switch.

- Double click the *Test limit switches* button. The window presented in Fig. 2.7 opens:




Fig. 2.7 Test limit switches window



Fig. 2.8 Switch detected window

Turn on manually one by one all switches related to  $x$ ,  $\theta$  and  $z$  directions. The  $z$ -axis switch turn on by your finger. The contactless switches turn on putting a coin close to the switches. When a switch is turned on one can hear a sound signal. If one turns on the  $x$ -axis limit switch then the window shown in Fig. 2.8 appears. It means that the switch operates properly. Close the window – click the *OK* button. When an arbitrary switch is not detected please check cables corresponding to the undetected switch.

Next, you can check if the cart, arm and payload move in the right direction and if the system stops at the desired limit position. The system is moved in the chosen direction until it reaches the zero position (at this point the switch limit must be active).

- Double click the *Go Home X-axis (T-angle direction and Z-axis)* button and observe the behaviour of the system. The window (Fig. 2.9) opens. You can interrupt the motion clicking the *OK* button.



Fig. 2.9 Go Home T-angle direction window

After performing tests along three directions the system is stopped at the zero position. The encoders of the  $x$ ,  $\theta$  and  $z$  directions are automatically reset to zero value.

If motion in a given direction is not observed check the cables and plugs related to this direction.

The next three steps perform the change of the system position from the initial position to the initial + 0.3 [m] or 0.3 rad position along a selected direction.

- Double click the *x-axis (  $\theta$  direction and z-axis) Movement* button. The window (see Fig. 2.10) opens where you can stop the motion by clicking the *OK* button.




Fig. 2.10 *x-axis movement* window

Click the *OK* button and the plot of the movements appears (see Fig. 2.11).




Fig. 2.11 Plot of the *x-axis* test movement

- In the next step double click the *Go To Centre* button. The system moves to the centre of the physical system workspace. The operational space boundaries are limited by the sizes of the laboratory set. They are fixed in the program. The window shown in Fig. 2.12 opens.




Fig. 2.12 *Go to Centre* window

Click the *OK* button the plot of the movement is displayed (Fig. 2.13)




Fig. 2.13 Plot of the movement to the center

Notice that the centre point has not been exactly reached. This is due to the open loop control mode. The control signal is turned off when the system exceeds the centre point.

The following steps are related to angle measurements.

- Double click the *Reset Angle Encoders* button. The window shown in Fig. 2.14 opens.




Fig. 2.14 *Reset Angle Encoders* window

Now you must set the load motionless and click *Yes*. The angle encoders are reset and the zero position is memorised by the system.

- To check if the angles are correct double clicks the *Check Angles* button. The window shown in Fig. 2.15 opens. Then manually move the load to a non-zero position and release it. Click the *OK* button. Observe the motion at the screen (Fig. 2.16)




Fig. 2.15 Start observation window




Fig. 2.16 Y vs. X angles trajectory

### 2.3. Stopping procedure

The system is equipped with the hardware stop pushbutton. It cuts off the transfer of control signals to the crane. The pushbutton does not terminate the real-time process running in the background. Therefore to stop the task you have to use *Simulation/Stop* from the pull-down menu in the model window.

### 3. Main Control Window

The user has a rapid access to all basic functions of the Tower Crane control system from *USB\_TCrane\_Main*. It includes tests, drivers, models and application examples.

In the *USB\_TCrane\_Main* window shown in Fig.2.4 the following groups may be distinguished

- Tools
- Device Driver and Simulation Model
- Demo Controllers
- Experiments

#### 3.1. Tools

The respective buttons in the TOOLS column perform the following tasks:

*Basic Tests* - checks the measurements and control.

*Go Home* – moves the crane to the zero position, resets the encoders and sets control signals to zero. This button is frequently used before starting an experiment. When the *Go Home* procedure is finished we can be sure that the values of all measured signals have been set to zero.

*Reset Angles* – resets the angle measuring encoders in a fixed position. If you stop the payload manually, perform the *Reset Angles* operation to be sure that the payload angles measured by the encoders show zeros.

*Go to Center* - moves the crane to the center of the crane operational space and switches off the control. Remember that the zero position of the crane is in the point where: the  $x$  position and the  $\theta$  angle are equal zero. Most experiments cannot be started from the zero point. *Go to Center* enables the crane quickly move to the center.

*Set Parameters* - enables the user to change the default values of *Rail limits*, *Base Address* and *z displacement*. The default value of the Base Address may cause a conflict with other devices installed in the computer. One has to be ensured that his computer configuration is free from address conflicts.

The user can also need to adjust the crane operational space dimensions to his requirements. Fig. 3.17 presents the window where such changes can be done. The user has to type numerical values into the editable text boxes.




Fig. 3.17 *Set Parameters* window



**All introduced modifications are written to the configuration file. Please be careful with introducing them. Writing the rail limit values exceeding the real rail limits may result in damage of the system elements.**

*Manual Setup* – The **Tower Crane Manual Setup** program gives access to the basic parameters of the laboratory 3-dimensional setup. The most important data transferred from the RT-DAC/USB board and the measurements of the crane signals as well as status signals and flags may be shown. Moreover the control signals of three DC drives may be set.

Double click the *Manual Setup* button and the screen shown in Fig. 3.18 opens.

The application contains five frames:

- The **RT-DAC/USB board** frame presents the main parameters of the USB board.
- The **Control** frame allows to change the control signals.
- The  $x$ ,  $\theta$  and  $z$  positions are given in the **X,  $\theta$  and Z positions** frame.
- The **X and Y angles** frame contains the angle measurements.
- The **Status and flags** frame displays state of the status signals and flag values.

All the data presented in the **Tower Crane Manual Setup** program are updated 20 times per second.






Fig. 3.18 The *Manual Setup* window

### RT-DAC/PCI / RT-DAC/USB board

The frame contains the parameters of the RT-DAC/PCI or RT-DAC/USB boards detected by the computer. With respect to the interface board applied to control the system the program operates with RT-DAC/PCI or RT-DAC/USB boards.

#### *No of detected boards*

The number of detected RT-DAC boards. If the number is equal to zero it means that the software has not detect any board. When more then one board is detected the *Board* list must be used to select the board that communicates with the program.

#### *Board*

The list applied to select the board currently used by the program. The list contains a single entry for each RT-DAC/PCI or RT-DAC/USB board installed in the computer. A new selection executed at the list automatically changes values of the remaining parameters within the frame. If more then one RT-DAC/USB board is detected the selection at the list must point to the board applied to control the tower crane system. Otherwise the program is not able to operate in a proper way.

#### *Bus number*

The number of the PCI bus where the current RT-DAC/PCI board is plugged-in. The parameter may be useful to distinguish boards, when more than one board is used and the computer system contains more than a single PCI bus. This parameter is only active if the RT-DAC/PCI boards are applied.

#### *Slot number*

The number of the PCI slot where the current RT-DAC/PCI board is plugged-in. The parameter may be useful to distinguish boards, when more than one board is used. This parameter is only active if the RT-DAC/PCI boards are applied.

#### *Base address*

The base address of the current RT-DAC/PCI board. The RT-DAC/PCI board occupies 256 bytes of the I/O address space of the microprocessor. The base address is equal to the beginning of the occupied I/O range. The I/O space is assigned to the board by the computer system and may be different among computers. This parameter is only active if the RT-DAC/PCI boards are applied. The base address is given in the decimal and hexadecimal forms.

#### *Logic version*

The number of the configuration logic of the on-board FPGA chip. A logic version corresponds to the configuration of the RT-DAC board defined by this logic and depends on the version of the tower crane model.

#### *I/O driver status*

The status of the driver that allows the access to the I/O address space of the microprocessor. The status has to be *OK* string. In other case the tower crane software HAS TO BE REINSTALLED.

## **Control**

The frame enables to set the control signals of three DC drives.

#### *X control, $\theta$ control, Z control*

The control signals of the *X*,  $\theta$  and *Z* DC drives may be set by entering a new value into the corresponding edit field or by dragging the corresponding slider. The control values may vary from  $-1.0$  to  $+1.0$ . The value of  $-1$ ,  $0.0$  and  $+1$  mean respectively: the maximum control in a given direction, zero control and the maximum control in the opposite direction to that defined by  $-1$ . If a new control value is entered in an edit field the corresponding slider changes respectively its position. If a slider is moved the value in the corresponding edit field changes as well.

#### *STOP*

The pushbutton is applied to switch off all the control signals. When pressed all the control values are set to zero.

#### *PWM prescaler*

The control signals are generated as PWM waves. The PWM prescaler sets the divider of the PWM reference signal. The frequency of the PWM controls is equal to:

$$F_{PWM} = 20000/1023/(1+PWMPrescalerr) [kHz]$$

This parameter sets the frequency of all PWM control signals. The parameter may vary from 0 to 63. It causes the changes of the frequency of the PWM control signals from 19.55kHz to 305Hz.

### ***X, $\theta$ and Z positions***

The frame presents data related to the  $X$ ,  $\theta$  and  $Z$  axis positions. The  $X$  position is the arm position. The  $\theta$  position relates to the angle position of the crane arm with the cart. The position of the load is denoted as  $Z$ . All position measurements are performed by the incremental encoders. There are the following four fields associated with each axis.

#### *Scale coefficient*

The value applied to calculate the position expressed in meters or radians. The value read from the encoder counter is multiplied by the corresponding scale coefficient to obtain the position in meters or radians.

#### *Position [bit]*

The value read from the corresponding encoder counter.

#### *Position [m]*

The position expressed in meters. The value read from the encoder counter is multiplied by the corresponding scale coefficient to obtain the position in meters.

#### *Position [rad]*

The position expressed in radians. The value read from the encoder counter is multiplied by the corresponding scale coefficient to obtain the position in radians.

#### *Reset*

The checkbox applied to reset the corresponding encoder counter. If the box is checked the related position is set to zero. The box has to be unchecked to allow position measurements.

As the incremental encoders are not able to detect an origin position the origin of the system has to be set by limit switches (see the description of the *Home autoreset flag*) or has to be set in a programming manner by a user. The **Reset** checkboxes are applied to set the origin position (zero position) in a programming way.

### **X and Y angles**

The frame presents data related to the  $X$  axis angles of the load. The angle measurements are performed by the incremental encoders. There are the following four fields associated with each angle.



***Angle  $X$  and angle  $Y$  denote the angle deviations in  $X$  and perpendicular to  $X$  direction, respectively. They are denoted by  $\alpha$  and  $\beta$  in the mathematical model (see Fig. 7.63).***

#### *Scale coefficient*

The value applied to calculate the angle expressed in radians. The value read from the encoder counter is multiplied by the corresponding scale coefficient to obtain the angle in radians.

*Angle [bit]*

The value read from the corresponding encoder counter.

*Angle [rad]*

The angle expressed in radians. The value read from the encoder counter is multiplied by the corresponding scale coefficient to obtain the position in radians.

*Reset*

The checkbox applied to reset the corresponding encoder counter. If the box is checked the related angle is set to zero. The box has to be unchecked to allow angle measurements.

As the incremental encoders are not able to detect an origin position the origin of the system has to be set in programming manner by a user. The angle reset checkboxes should be checked when the load remains motionless in the downright position.

## **Status and flags**

The frame presents status data and flags related to the  $X$ ,  $\theta$  and  $Z$  directions. There are seven fields associated with each directions.

*Position limit [64\*bit]*

The board is able to automatically turn off the control signal if the tower crane system is going outside the operating range. The field defines the maximum value of the corresponding encoder determining the maximum accessible position. If the encoder reaches this position the board is able to stop the generation of the DC control signals moving the system outside the operating range.

*Position limit [m] or [rad]*

The field defines the maximum accessible position expressed in meters or radians. If the system reaches this position the board is able to stop the generation of the DC control signals moving the system outside the operating range. The maximum position expressed in meters or angle position expressed in radians is obtained as the result of the multiplication of 64, maximum position expressed in bits and the corresponding scale coefficient.

*Max position limit flag*

If the checkbox is selected the board turns off the control when the system is going to move outside the operating range. It is recommended to keep this checkboxes selected all the time.

*Home position limit switch*

The boxes that present the state of the home limit switches. If a limit switch is pressed the corresponding box is selected.

*Home autoreset flag*

The flag that causes automatic reset of the encoder counter when the corresponding home limit switch is pressed. If the checkbox is unchecked (the flag is inactive) the state of the limit switch does not influence the state of the encoder counter.

#### *Therm state*

The signal that presents the state of the thermal flag of the power interface. The system contains three power amplifiers for the DC drives. If the power amplifier is overheated the corresponding box is checked.

#### *Therm flag*

The flag that causes to turn off the control if the power amplifier is overheated. If the flag is unchecked, the power amplifier is overheated and the temperature increases the amplifier itself turns off the control signal. It is recommended to keep this checkboxes selected all the time.

Click the *Go to Center* button and open next *Manual Setup* button. You will see the changes of measured positions. Note that the X position,  $\theta$  angle position and Z position have changed their values and become the center positions in the crane workspace.

### 3.2. Drivers

The main driver is located in the *RTWT Device Driver*. The driver is a software go-between for the real crane MATLAB environment and the RT-DAC/USB board. This driver is dedicated to the control and measurement signals. Click the *Tower Crane Device Drivers* button and the driver window will be opened (Fig. 3.19)




Fig. 3.19 Tower Crane Device Drivers

The driver has three PWM inputs (DC motor controls) for the X and Z-axes and T- *angle direction*. There are 10 outputs of the driver: X position, T position, Z position, two angles (see section 8) and additionally three safety switches. According to a pre-programmed logic the internal XILINX program of the RT-DAC/USB board can use the switches to stop the DC motors.

When one wants to build his own application one can copy this driver to a new model.



**Do not make any changes inside the original driver. They can be made only inside its copy!!!**

### 3.3. Simulation Model

The simulation model of the crane is located under this button. Its signal environment is identical as the model given in the *RTWT Device Drivers* except the lack of the safety switches (see Fig. 3.20). These switches are not used in the simulation mode. The simulation model is used for many purposes: identification, controllers design, etc.




Fig. 3.20 Tower Crane Simulation Model

In the mask shown in Fig. 3.21 one can introduce initial values for the model state variables. Only the model with a constant value of the rope length is accessible, so the appropriate check box (Fig. 3.21) has to be checked. (see section 7.1 for details).




Fig. 3.21 Mask of the simulation model

The simulation model is running in the normal simulation mode.



**Set solver options to: *Fixed step* and *Fixed-step size* equal to 0.001 atmost. A greater value results in errors.**

The C source code of the *modeltc.c* file is attached to the *DevDriv* directory.

### 3.4. Demo Controllers

A number of control algorithms are given. These demos can be used to familiarise the user with the crane system operation and enable to create the user-defined control systems. The examples must be rebuilt before using. Due to similarity of the examples we focus our attention only on one of them.

Click the *Relay* button the model appears (see Fig. 3.22):




Fig. 3.22 Control system with the relay controllers.

Notice that the model looks like a typical Simulink model. The device driver is applied in the same way as other blocks from the Simulink library. The only difference is that the model is used by RT-CON to create the executable library, which runs in the real-time mode.

The goal of the model is the relay control in the x-axis only. Therefore *GainY* and *GainZ* are set to zero.

- Look at the mask of the *Relay* block connected to the X PWM input (Fig. 3.23).  
Note that the control generated by the controller has two values:  $+0.5$  and  $-0.5$ . The On/Off limits are 0.2 and 0.6. This means that the crane will move between these limits with the speed corresponding to the control value equal to 0.5.
- To choose the starting point inside the  $[-0.5, 0.5]$  range go to the *USB\_TCrane\_Main* window and click the *Go to Center* button.
- Then, choose the *Tools* pull-down menus in the Simulink model window. The pop-up menus provide a choice between predefined items. Choose the *RTW Build* item. A successful compilation and linking process is finished with the following message:






Fig. 3.23 X-PWM controller parameters

*Model TowerCrane\_Relay.rtd successfully created*

*### Successful completion of Real-Time Workshop build procedure for model: TowerCrane\_Relay*

If any error occurs then the message corresponding to the error is displayed in the MATLAB command window. Next, click the *Tools/External Mode Control Panel* item and next click the *Signal Triggering* button. The window shown in Fig. 3.24 opens.

- Select *XT Scope* , set *Source* as the manual option, mark the *Arm when connect to Target* option and close the window.
- Return to the model window and click the *Simulation/Connect to Target* option. Next, click the *Simulation/Start real-time code* item.




Fig. 3.24 External Signal & Triggering window

- Observe the plots in the scope and click *Stop Simulation* after some time. Results of the example are shown in Fig. 3.25.




Fig. 3.25 Results of the relay controller demo experiment

The  $X$  position starts from 0.38 and changes between 0.4 and 0.1. The control (magenta line) is a square wave in the range  $[-0.5, 0.5]$ . The control is switched when the limit values are reached by the  $X$  position. Note that the  $X$  angle in the form of a sinusoidal curve is modulated by the control interacting with friction.

#### 4. Your first real-time control experiment

We recommend two experiments. In the first one, the control loop only in the  $x$  direction is defined. In this case stabilisation of the angle of the payload is neglected. In the second experiment the stabilisation of the payload angle is added.

We begin from a simple real-time control experiment. A PID controller for the  $x$  position of the cart is built. The *USB\_TCrane\_first* model is given in Fig. 4.26. To invoke it, click *Model for control experiment* button in the *The USB\_TCrane\_Main* window. In this case, the active control corresponding to the  $x$  cart position is all you need. The Simulink blocks included in the control are drawn with dropped shadows to distinguish active control loops from disabled loops. In our first experiment we use the *X position of the cart* PID controller with its P control part activated only.




Fig. 4.26 Two PID controllers applied in a real time experiment. The first step – PID of *X position of the cart* is active. The second step – both controllers are active

We define a source of a desired *x* position signal as the *X reference/Generator* from the Simulink library (see Fig. 4.27)




Fig. 4.27 Signal Generator becomes the square wave signal source

As usual, at the beginning the *GO HOME* and *GO TO CENTER* actions are performed. The proper experiment starts in the middle of the  $x, y$  rails. Therefore the constant block *shift* is required.

Next, we define the signals to be used. These are *X-reference* – the desired  $x$  position value, *X-position*, *X-control* and *X-angle*. These signals, among others, are connected to the *Scope* block.

The properties of this block are defined below (see Fig. 4.28). This window opens after the selection of the *Scope/Properties* tab. Mark the *Save data to workspace* checkbox, define the *Variable name* as *EX1* and the data format as structure. This means the collected data within 30 seconds time range are saved to the workspace in the structure *EX1*. The sampling period set in the *Simulation Parameters* window (see Fig. 4.28) is equal to 0.01, *Sampling Decimation* is set to 10. Therefore, the size of *EX1* is equal to:  $30\text{ s} / (0.01\text{ s} \times 10) + 1 = 301$




Fig. 4.28 Setting of the *Scope* block

Next, return to *USB\_TCrane\_first* and select the *Simulation/Parameters* item. In the *Solver* tab select *Fixed Step* and set *Stop time* equal to 30. The *Real-time Workshop* options must be defined as in Fig. 4.29.




Fig. 4.29 The *Configuration Parameters* window

Next connection with external mode and transport layer has to be set. Interface page, responsible for connection with external target system, is shown in Fig. 4.30. To fulfil it properly a user must select *External mode* in *Interface* edit window. After that the user have to set *RT-CON tcpip* in the *Transport layer* edit window.




Fig. 4.30 External mode settings

Set the PID controllers. We set the *Proportional* part of the *X position of the cart* PID controller as indicated by the arrow in Fig. 4.31. The *X angle of the payload* PID controller remains inactive. The other controlling loops are disabled due to the *Gain* blocks set to zero (see Fig. 4.26).




Fig. 4.31 Setting of the PID controllers

Mark *Simulation/External* item in the *TowerCrane\_first* model window (see Fig. 4.32).




Fig. 4.32 External control mode

Next, invoke the *Tools/External mode control panel* item. The *External Mode Control Panel* window opens (see Fig. 4.33).




Fig. 4.33 Setting data acquisition options

By clicking on the *Signal & Triggering* button invoke the window shown in Fig. 4.34. In this window we define a triggering mode for marked blocks. In our case only one block exists – *XT Scope*. We mark *XT Scope*, set *Source* as the *manual* option, mark *Arm when connect to target* and close the window.

Now we can build the real-time model. To do it select the *Simulation/Parameters* item and then the *Real-Time Workshop* option in the model window, and click the *Build* item. Successful compilation and linking process should be finished with the following message:

```
### Created RT-CON executable: USB_TowerCrane_first.dll
```

```
### Successful completion of Real-Time Workshop build procedure for model:
USB_TowerCrane_first
```




Fig. 4.34 Setting triggering of signals

## 4.1. Real-time experiment

Having prepared the controller model you can start the real-time experiment. Two actions *GO HOME* and *GO TO CENTER* must be performed at the beginning in the *USB\_TCrane* window. The crane is ready for the experiment. The cart is in the middle of the cylindrical  $x$ ,  $\theta$  plane, the payload is hanging down in its rest position. Open the *Scope* figure clicking on the *Scope* block.

Now return to the model window and click the *Simulation/Connect to Target* item. Next, click the *Simulation/Start real-time code* item. It activates the experiment lasting 30 seconds. Observe the cart motion in the  $x$  direction. The cart follows the desired square wave signal controlled by the P regulator. The payload oscillates freely, being uncontrolled. After 30 seconds the experiment stops. The history of the *EX1* variable is visible in the *Scope* (see Fig. 4.35). Notice the harmonic (uncontrolled) angle signal of the payload, the square wave generated by *Signal Generator* followed by the cart  $x$  position signal. The static error is due to the inadequate P control action. The control has the highest magnitude among other signals visible in the figure. When an abrupt change in the wave signal occurs then it results in the saturation of the control signal.




Fig. 4.35 Data visible in the scope during the experiment

## 4.2. Data processing

The results are saved to the workspace as a structure variable *EX1*. If you write the variable name in the MATLAB command window then you obtain the answer

```
EX =
    time: [301x1 double]
   signals: [1x1 struct]
  blockName: 'TCrane_pid_all/Scope'
```

This data can be plotted in many ways. For example use the following command

```
>> plot(EX1.time,EX1.signals.values(:,1:4))
```


You can repeat the experiment several times using different P parameter settings and including another P controller for the *x* angle. The parameters of the controllers for successive five experiments are given in Table 4-1.

Number of experiment	P of the cart position	P of the payload angle	Figure
1	2.5	—	Fig. 4.36
2	5	—	Fig. 4.36
3	2.5	1	Fig. 4.37
4	5	2	Fig. 4.37
5	5	4	Fig. 4.37


Table 4-1 Parameters for the experiments



The results of five experiments are presented in Fig. 4.36 and Fig. 4.37.




P set to 2.5




P set to 5


Fig. 4.36 Only one controller is active – desired  $x$  position of the cart is tracked



P of the cart position set to 2.5  
P of the payload angle set to 1



P of the cart position set to 5  
P of the payload angle set to 2



P of the cart position set to 5  
P of the payload angle set to 4

Fig. 4.37 Two controllers – desired  $x$  position of the cart and  $\mathcal{X}$  angle of the payload are tracked

The left-hand side of Fig. 4.36 shows that the P value set to 2.5 is too small for a proper  $x$  position tracking. The static error of  $x$  position is large. A higher gain P equal to 5 (the right-hand side of Fig. 4.36) reduces the static error but results in the saturation of control.

In Fig. 4.37 one can see similar results obtained for two active controllers. We focus our attention on the  $x$  angle control. The trade-off between two acting controllers is well visible in the upper left-hand picture of Fig. 4.37. One control signal serves for two control purposes: follow the desired value of the cart  $x$  position and simultaneously stabilises the payload in its hanging down position.

## 5. Simulation

### 5.1. Simulation control experiments

We can execute the experiment similar to experiments from the previous section in a purely simulated form.

We invoke the *USB\_TCrane\_first\_model* window. Notice two differences. The *TCrane* real-time driver block has been replaced by the *Tower Crane model* simulation model block. The *External* mode of operation has been replaced by the *Normal* mode of operation (see Fig. 5.38).




Fig. 5.38 Two PID controllers applied in a simulated experiment. The first step – PID of  $X$  position of the cart is active.

The interior of *Tower Crane model* includes the complete non-linear model described in section 7. After clicking on *Tower Crane model* the mask given in Fig. 5.39 opens. You can set the initial values of eight variables. The length of the lift line is constant. In our example all initial conditions are set to zero. The  $Z$  position is set to 0.6 and the last variable  $t$  is set to -1. Setting  $t$  to -1 denotes that the source of time is the *RT-CON* clock.




Fig. 5.39 Mask of *Tower Crane model*

If you look under the *Tower Crane model* mask you can see its interior (Fig. 5.40). *modeltcsimple* is the executable dll file. The scale factors are set to 1. These parameters relate to identification parameters like friction and tension of the belts.




Fig. 5.40 The interior of *Tower Crane\_model*

In opposite to the real-time model it is not necessary to rebuild the model before running. Open *Simulation parameters* from the menu bar. Notice the solver options: *Fixed step* and *Fixed-step size* set to 0.001 in the editable text box. You can start a simulation from the menu bar. When the simulation is running you can watch the results in the scope (see Fig. 5.41). Similarly as in the real time experiment, the data are saved in the scope in the *EX* structural variable.

Plotting four curves as functions of time (Fig. 5.41) is produced by:

```
>> plot(EX.time, EX.signals.values)
```

We perform simulations related to the previous experiments. The parameters of the *xw pos* PID regulator are set to:  $P=3$ ,  $I=0.3$ . The results are visible in Fig. 5.41. They are similar to the results of the real-time experiments. The model reflects characteristic features of the laboratory crane. The model compatibility to the real crane strongly depends on such parameters as the payload lift-line length, the static cart friction, the belt tension, etc. These and other parameters are written into the C source code of the *modeltowecrane.c* file attached in the *DevDriv* directory. If you wish to modify this file please make a copy and introduce the new parameters in the copy. Remember to produce an executable .dll file afterwards.




Fig. 5.41 Simulation data visible in the scope

Let us turn on the controller PID *X ag* - responsible for dumping the *X* angle oscillation. To do this, we change the value of the gain vector ( “C” block ) from the  $[0 \ 1]$  value to the  $[1 \ 0.3]$ , see Fig. 5.42.




Fig. 5.42 New parameters of the “C” block

The calculations of the dumping oscillations controller  $X_{ag}$  are taken now with the weight equal 1 while the  $xw_{pos}$  PID controller acts with the weight 0.3. This causes, the the control algorithm mainly dumps the oscillation. The result of the experiment with the working  $X_{ag}$  controller presents Fig. 5.43.



Fig. 5.43 Dumped oscillations of the X-angle

One may observe, that the oscillations are minimized by the controller, but a desired position of the cart (square wave) is tracking in a worse way then in the previous experiment.

## 5.2. PID simulation control of load position

In our experiment the cart is following the desired trajectory similar to Lissajous curves. We invoke the *TCrane\_impres* model from MATLAB Command Window (see Fig. 5.44). We put the cart into motion in two directions. The desired cart positions are generated as two sinusoidal signals. There are two generators identical in amplitudes equal to 0.2 m and different in frequencies. The *X* motion frequency is 0.05 rad/s (Fig. 5.45) and the *T* angle motion frequency is 0.1 rad/s.




Fig. 5.44 The controller built for the cart to follow a Lissajous curve

We perform the experiment twice, for first time without two angle controllers (*X ag* and *Y ag*). For the second time the angle controllers are included. The weights of each controller are set in the *C* and *C1* Simulink blocks. For the first experiment we set the gains to the  $[0 \ 1]$  values (see Fig. 5.45 and Fig. 5.46). For the second experiment we set the gains of the *C* and *C1* block to the  $[0.3 \ 1]$  values.




Fig. 5.45 Gain value of the *C* block




Fig. 5.46 Gain value of the *C1* block

The cart motion is shown in Fig. 5.47. The thick line represents the cart position in the *XY* plane (there is no movement in the *Z* direction). The respective controls and the payload angles are shown in Fig. 5.48 and Fig. 5.49.




Fig. 5.47 The desired cart positions – thin line and the cart positions – thick line (in meters)




Fig. 5.48 The controls (normalised units) and the payload angles [rad] on the  $XY$  plane; the angles without control




Fig. 5.49 The controls (normalised units) and the payload angles [rad] on the  $XY$  plane; the angles with control

Notice that the control curve in the  $XT$  plane for the case of uncontrolled deviations of the payload is smoother than that for the controlled deviations case. The payload deviations are zoomed in Fig. 5.50 and Fig. 5.51 (please notice the scales of deviations).




Fig. 5.50 Unstabilised deviations of the payload in the  $X Y$  plane (in radians)




Fig. 5.51 Stabilised deviations of the payload in the  $X Y$  plane (in radians)

The controls in the  $X$  and  $T$  directions are shown in Fig. 5.52 and Fig. 5.53. Two cases of control: with and without stabilisation of angles are compared. The magnifications of the chosen areas are present in Fig. 5.54 and Fig. 5.55




Fig. 5.52 The controls (normalised units) in the  $X$  direction vs. time (seconds)




Fig. 5.53 The controls (normalised units) in the  $T$  angle direction vs. time (seconds)

Notice that the controls related to the case of angles stabilisation share their actions between two tasks: tracking a desired position of the cart and stabilising the payload in its down position.




Fig. 5.54 The controls (normalised units) in the  $X$  direction vs. time (seconds)




Fig. 5.55 The controls (normalised units) in the  $T$  angle direction vs. time (seconds)



## 6. Prototyping your own controller in the real-time environment

In this section the process of building your own control system is described. The *RTW* toolbox of MathWorks and RT-CON toolbox of INTECO are used. In this section we give indications how to proceed in the real-time environment.

Before starting, test your MATLAB configuration and compiler installation by building and running an example of a real-time application. TCrane toolbox includes the real-time model of PC speaker.

The model `pc_speaker_xxx.mdl` does not have any I/O blocks so that you can run this model regardless of the I/O boards in your computer. Running this model will test the installation by running Real-Time Workshop, and your third-party C compiler.

In the MATLAB command window, type

**pc\_speaker\_vc**      (if you are using Visual C++ compiler)  
or  
**pc\_speaker\_watc**      (if you are using OpenWatcom 1.3 compiler)

Build and run this real-time model as follows:

- From the *Tools* menu, point to *Real-Time Workshop*, and then click *Build Model*.

The MATLAB command window displays the following messages.



```
#### Starting Real-Time Workshop build procedure for model: pc_speaker_vc
#### Generating code into build directory: C:\apps\MATLAB\R2006a\work\
pc_speaker_vc_RTCON
#### Invoking Target Language Compiler on pc_speaker_vc.rtw. . .
.
. . .
#### Created RT-CON executable: pc_speaker_vc.dll
#### Successful completion of Real-Time Workshop build procedure for model:
pc_speaker_vc
```

- From the *Simulation* menu, click *External*, and then click *Connect to target*.

The MATLAB command window displays the following message.  
Model `pc_speaker_vc` loaded

- From *Simulation* menu, click *Start real-time code*.

The *Dual* scope window displays the output signals and you hear the speaker sound.

Only correct configuration of the MATLAB, Simulink, RTW and C compiler guaranties the proper operation of the Tower Crane system.

To build the system that operates in the real-time mode the user has to:

- create a Simulink model of the control system which consists of the *Tower Crane Driver* and other blocks chosen from the Simulink library,
- build the executable file under RTW and RT-CON (see the *Tools* menu in Fig. 5.56)
- start the real-time code using menu of the model window.




Fig. 5.56 Creating the executable file using RTW and RT-CON

## 6.1. Creating a model

The simplest way to create a Simulink model of the control system is to use one of the models included in the *USB\_TCrane\_Main* window as a template. For example, click on the *Relay* button and save it as *MySystem.mdl* name. The *MySystem* Simulink model is shown in Fig. 5.57.




Fig. 5.57 The *MySystem* Simulink model

Now, you can modify the model. You have got absolute freedom to develop your own controller. Remember to leave the *Tower Crane Device* driver model.

Though it is not obligatory, we recommend you to leave the multiplexer with the scope and the control saturation blocks. You need a scope to watch how the system runs. You also need the saturation blocks to constraint the controls to match the maximal PWM signals sent to the DC motors. The saturation blocks are built in the *Tower Crane* driver block. They limit currents to the DC motors for safety reasons. However they are not visible for the user who can be surprised that the controls saturate. Other blocks that remain in the window are not necessary for our new project.

Creating your own model on the basis of an old example ensures that all-internal options of the model are set properly. These options are required to proceed with compiling and linking in a proper way. To put the *Tower Crane Device Driver* into the real-time code a special make-file is required. This file is included to the *Tower Crane* software.

You can apply most of the blocks from the Simulink library. However, some of them cannot be used (see MathWorks references manual).

The scope block properties are important for appropriate data acquisition and watching how the system runs.

The *Scope* block properties are defined in the Scope property window (see Fig. 5.58). This window opens after the selection of the *Scope/Properties* tab. You can gather measurement data to the *Matlab Workspace* marking the *Save data to workspace* checkbox. The data is placed under *Variable name*. The variable format can be set as *structure* or *matrix*. The default *Sampling Decimation* parameter value is set to 1. This means that each measured point is plotted and saved. Often we choose the *Decimation* parameter value equal to 5 or 10. This is a good choice to get enough points to describe the signal behaviour and to save the computer memory.




Fig. 5.58 Setting the parameters of the *Scope* block

When the Simulink model is ready, click the *Tools/External Mode Control Panel* option and next click the *Signal Triggering* button. The window presented in Fig. 5.59 opens. Select *XT Scope*, set *Source* as manual, set *Duration* equal to the number of samples you intend to collect and close the window.




Fig. 5.59 External Signal & Triggering window

## 6.2. Code generation and the build process

Once a model of the system has been designed the code for the real-time mode can be generated, compiled, linked and downloaded into the processor.

The code is generated by the use of Target Language Compiler (TLC) (see description of the Simulink Target Language). The make-file is used to build and download object files to the target hardware automatically.

First, you have to specify the simulation parameters of your Simulink model in the *Configuration parameters* dialog box. The RTW page appears when you select the *Real-TimeWorkshop* option (Fig. 5.60). The *RTW* page enables to set the real-time build options and then to start the building process of the *RTW.DLL* executable file.




Fig. 5.60 RTW option of the *Configuration parameters* dialog box

The system target file name is *towercrane.tlc* (for M7 R2006a/b R2007a/b, R2008a/b). It manages the code generation process. The *towercrane\_vc.tmf* template makefile is responsible for C code generation using the open Watcom compiler.

The *Solver* page appears when you select the *Solver* option (Fig. 5.61). The *Solver* page enables to set the simulation parameters. Several parameters and options are available in the window. The *Fixed-step size* editable text box is set to 0.01 (this is the sampling period in seconds).



**The *Fixed-step* solver is obligatory for real-time applications. If you use an arbitrary block from the discrete Simulink library or a block from the driver library remember that different sampling periods must have a common divider.**




Fig. 5.61 Solver option of the Configuration parameters page

The *Start time* has to be set to 0. The solver has to be selected. In our example the fourth-order integration method – *ode4(Runge-Kutta)* is chosen.




Fig. 5.62 Interface option of the Configuration parameters page

If all parameters are set properly you can start the DLL executable building process. For this purpose press the *Build* push button on the RTW page (Fig. 5.60) or click Ctrl-B at

keyboard being in the model window.. Successful compilation and linking processes generate the following message:

*Model MyModel.rtd successfully created*

*### Successful completion of Real-Time Workshop build procedure for model: MyModel*

Otherwise, an error message is displayed in the MATLAB command window.



**Before starting the experiment set the initial position of the cart, the arm and the payload in a safe zone. The *Go Home* and *Go To Center* buttons are applied to fulfil these tasks.**



## 7. Mathematical model of the Tower Crane

### 7.1. Equations

The schematic diagram related to the tower crane mathematical model is shown in Fig. 7.63.




Fig. 7.63 Tower Crane model coordinates

There are five measured quantities:

- $x_w$  denotes the distance of the cart along the arm from the center of the construction frame;
- $\theta$  denotes the angular position of the arm
- $L$  denotes the length of the lift-line;
- $\alpha$  denotes the angle between the  $z$  axis and the projection of the lift-line onto  $xz$  plane;
- $\beta$  denotes the angle between the projection of the lift-line onto the  $xz$  plane and the lift-line;
- $x_c, y_c, z_c$  define coordinates of the payload.

An important element in the construction of the mathematical model is the appropriate choice of the coordinates. The center point ( $x=0, y=0, z=0$ ) of the Cartesian system is in the center of the crane tower at the level equal to the suspension point of the load. The position of the payload is described by two angles,  $\alpha$  and  $\beta$ , shown in Fig. 7.63. The Cartesian system rotates according to the tower movement, described by the  $\theta$  angle.

It is assumed that the lift-line is permanently stretched. The position of the payload according to Fig. 7.63 is described by the formulas:

$$x_c = x_w - L \cos \beta \sin \alpha \quad (1)$$

$$y_c = L \sin \beta \quad (2)$$

$$z_c = -L \cos \alpha \cos \beta \quad (3)$$

The first derivatives of the equations (1-3) and the kinetic and potential energy of the payload are calculated. Based on the Lagrange equations, Amjed A. Al-Mousa<sup>1</sup> gives the following equations that describes the dynamics of the tower crane.

$$\begin{aligned}
\dot{x}_1 &= x_5 \\
\dot{x}_2 &= x_6 \\
\dot{x}_3 &= x_7 \\
\dot{x}_4 &= x_8 \\
\dot{x}_5 &= -\frac{1}{2L} \left( 2g \cos x_2 \sin x_1 + 4x_7 x_8 \cos x_1 - L x_8^2 \sin(2x_1) \cos^2 x_2 - 2x_3 x_8^2 \sin x_1 \sin x_2 \right. \\
&\quad \left. - 4L x_6 x_8 \cos x_2 \cos^2 x_1 + L x_6^2 \sin(2x_1) + 2 \sin x_1 \sin x_2 u_1 + 2x_3 \cos x_1 u_2 \right. \\
&\quad \left. - 2L \sin x_2 u_2 \right) \\
\dot{x}_6 &= -\frac{1}{L \cos x_1} \left( g \sin x_2 + x_3 x_8^2 \cos x_2 - L x_8^2 \sin x_2 \cos x_1 \cos x_2 + 2L x_5 x_8 \cos x_1 \cos x_2 \right. \\
&\quad \left. - 2L x_5 x_6 \sin x_1 - \cos x_2 u_1 + L \sin x_1 \cos x_2 u_2 \right) \\
\dot{x}_7 &= u_1 \\
\dot{x}_8 &= u_2
\end{aligned}$$

where:  $x_1 = \beta$  ,  $x_2 = \alpha$  ,  $x_3 = x_w$  ,  $x_4 = \theta$  ,  $x_5 = \dot{\beta}$  ,  $x_6 = \dot{\alpha}$  ,  $x_7 = \dot{x}_w$  ,  $x_8 = \dot{\theta}$  ,  
 $u_1 = \ddot{x}_w$  ,  $u_2 = \ddot{\theta}$  .

It is assumed that the length of the lift line  $L$  is constant.

## 7.2. Comparison between the mathematical model and laboratory system.

Figures 6.64 and 6.65 presents variables and directions used in the model and system




Fig. 7.64. Variables and axis assumed at the mathematical model




Fig. 7.65. Variables and axis adopted at the mechatronic laboratory system

The basic differences between the model and system are as follows:

- the control signal – it is the acceleration inside the model it is the velocity inside the system,
- the angles of the payload – see the completely different notification in the above figures.
- the tower rotate clockwise in the model and counterclockwise in the system.

<sup>1</sup> Amjed A. Al-Mousa, *Control of Rotary Cranes Using Fuzzy Logic and Time-Delayed Position Feedback Control*, Virginia Polytechnic Institute and State University, 27-th November 2000

The Tower Crane system is supplied with a Simulink block *Tower Crane model* that is based on the mathematical model see Fig. 7.66.




Fig. 7.66. Simulink representation of the mathematical model

The model is equipped in two control inputs:  $X$  control and  $T$  control that correspond to the appropriate accelerations. The third input feeds a constant as the lift line length. The state variables are collected as the outputs.

Double clicking the *Tower Crane model* block to introduce the initial values of the states Fig. 7.67.




Fig. 7.67. Simulink window to introduce the initial values of the states

To adopt the mathematical model to the system a number of modifications of the Simulink model have been introduced (see Fig. 7.68).




Fig. 7.68. The modified Simulink window to introduce the initial values of the states

To assure a control signal compatibility (velocity – acceleration) two Simulink block *Derrivative* are added and two gains: *Scale rate x* and *Scale rate T*.

A conversion between the angles  $\alpha$ ,  $\beta$  and  $X, Y$  describes two formulas:  

$$X = \arcsin(-\cos \beta \cdot \sin \alpha)$$
  

$$Y = \arcsin\left(-\frac{\sin \beta}{\cos X}\right)$$

The clockwise rotation in the system and counterclockwise rotation in the model have been equalized. The gain = - 1 in the  $n2$  block was introduced and the Scale rate T factor was multiplied by - 1 as well.

## 8. Description of the Tower Crane class properties

The *towercrane* is a MATLAB class, which gives the access to all the features of the RT-DAC/USB board equipped with the logic for the Tower Crane model. The RT-DAC/USB board is an interface between the control software executed by a PC computer and the power-interface electronic of the Tower Crane model. The logic on the board contains the following blocks:

- incremental encoder registers – five 16-bit registers to measure the position of the incremental encoders. There are five encoders measuring five state quantities: two cart positions at the cylindrical plane, the lift-line-length and two deviation angles of the payload;
- incremental encoder resets logic. The incremental encoders are able to generate different output waves when the encoder rotates clockwise and the counter clockwise. The encoders are not able to detect the reference (“zero”) position. To determine the “zero” position the incremental encoder registers can be set to zero from the computer program or an encoder register is reset when the corresponding limit switch to the encoder is reached;
- PWM generation block – generates three sets of signals. Each set contains the PWM output signal, the direction signal and the brake signal. The PWM prescaler determines the frequency of all the PWM waves. The PWM block logic can prevent the cart from motion outside the rail limits and the lift-line angles from lying outside the operating range. The operating ranges are detected twofold – by the limit switches and by three limit registers;
- power interface thermal flags – when the temperature of the power interface for the DC motors is too high the thermal flags can be used to disable the operation of the corresponding overheated DC motor.

All the parameters and measured variables from the RT-DAC/USB board are accessible by appropriate methods of the *towercrane* class.

The object of the *towercrane* class is created by the command:

```
object_name = towercrane;
```

The *get* method is called to read a value of the property of the object:

```
property_value = get(object_name, 'property_name');
```

The *set* method is called to set new value of the given property:

```
set(object_name, 'property_name', new_property_value);
```

The *display* method is applied to display the property values when the *object\_name* is entered in the MATLAB command window.

This section describes all the properties of the *towercrane* class. The description consists of the following fields:

Purpose	Provides short description of the property
Synopsis	Shows the format of the method calls
Description	Describes what the property does and the restrictions of is subjected to
Arguments	Describes arguments of the set method
See	Refers to other related properties
Examples	Provides examples how the property can be used

## 8.1. BaseAddress

**Purpose:** Read the base address of the RT-DAC/USB board.

**Synopsis:** *BaseAddress = get( tcrane, 'BaseAddress' );*

**Description:** The base address of RT-DAC/USB board is determined by OS. Each *towercrane* object has to know the base address of the board. When the *towercrane* object is created the base address is detected automatically. The detection procedure detects the base address of the first RT-DAC/USB board plugged into the USB slots.

**Example:** Create the *towercrane* object:

```
tcrane = towercrane;
```

Display its properties by typing the command:

```
tcrane
```

```
>>Type:          towercrane Object
>>BaseAddress:   528
>>Bitstream ver.: x33
>>Encoder:       [ 65479  7661  20032  65533  65534 ][bit]
>>               [ -0.0022207[m] 0.29847[m] 0.38411[m] -0.004602[rad] -0.003068[rad] ]
>>Z displacement: 0.32[m]
>>PWM:           [ -0.062561  0.031281      -1 ]
>>PWMPrescaler:  60
>>RailLimit:     [ 361  381  815 ]*64[bit] <--> [23104  24384  52160 ][bit]
>>               [ 0.90013    0.95    1.0002 ][m]
>>RailLimitFlag: [ 1  1  1 ]
>>RailLimitSwitch: [ 0  1  1 ]
>>ResetSwitchFlag: [ 0  0  0 ]
>>Therm:         [ 1  1  1 ]
>>ThermFlag:     [ 1  1  1 ]
>>Time:          1.041 [sec]
```

Read the base address:

```
BA = get( tcrane, 'BaseAddress' );
```

## 8.2. BitstreamVersion

**Purpose:** Read the version of the logic design for the RT-DAC/USB board.

**Synopsis:** *Version = get( tcrane, 'BitstreamVersion' );*

**Description:** This property determines the version of the logic design of the RT-DAC/USB board. The Tower Crane models may vary and the detection of the logic design version makes it possible to check if the logic design is compatible with the physical model.

### 8.3. Encoder

**Purpose:** Read the incremental encoder registers.

**Synopsis:** `enc = get( tcrane, 'Encoder' );`

**Description:** The property returns five digits. The first two measure the position of the cart – linear and angle. The third digit is used to measure the length of the lift-line and the last two measure the angles of the lift-line. The returned values can vary from 0 to 65535 (16-bit counters). When a register is reset the value is set to zero. When a rail limit flag is set it disables the movement outside the defined working range (rail limit). When a reset switch flag is set the encoder register is reset automatically when the appropriate switch is pressed.

The incremental encoders generate 4096 or 2048 pulses per rotation. The values of the *Encoder* property should be converted into physical units.

**See:** *ResetEncoder, RailLimit, RailLimitFlag, ResetSwitchFlag*

### 8.4. PWM

**Purpose:** Set the parameters of the PWM waves.

**Synopsis:** `PWM = get( tcrane, 'PWM' );`  
`set( tcrane, 'PWM', NewPWM );`

**Description:** The property determines the duty cycle and direction of the PWM waves for three DC motors. The first two DC motors control the position of the cart and the last motor controls the length of the lift-line. The *PWM* and *NewPWM* variables are 1x3 vectors. Each element of these vectors determines the parameters of the PWM wave for one DC motor. The values of the elements of these vectors can vary from -1.0 to 1.0. The value -1.0 means the maximum control in one direction, the value 0.0 means zero control and the value 1.0 means the maximum control in the opposite direction.

The PWM wave is not generated if:

- a rail limit flag is set and the cart or lift-line are going to operate outside the working range,
- a therm flag is set and the power amplifier is overheated.

**See:** *RailLimit, RailLimitFlag, Therm, ThermFlag*

**Example:** `set( tcrane, 'PWM', [ -0.3 0.0 1.0 ] );`

### 8.5. PWMPrescaler

**Purpose:** Determine the frequency of the PWM waves.

**Synopsis:** `Prescaler = get( tcrane, 'PWMPrescaler' );`  
`set( tcrane, 'PWMPrescaler', NewPrescaler );`

**Description:** The prescaler value can vary from 0 to 63. The 0 value generates the maximum PWM frequency. The value 63 generates the minimum frequency.

**See:** *PWM*

## 8.6. ResetEncoder

**Purpose:** Reset the encoder counters.

**Synopsis:** *set( tcrane, 'ResetEncoder', ResetFlags );*

**Description:** The property is used to reset the encoder registers. The *ResetFlags* is a 1x5 vector. Each element of this vector is responsible for one encoder register. If the element is equal to 1 the appropriate register is set to zero. If the element is equal to 0 the appropriate register remains unchanged.

**See:** *Encoder*

**Example:** To reset the first and fourth encoder registers execute the command:  
*set( tcrane, 'ResetEncoder', [ 1 0 0 1 0 ] );*

## 8.7. RailLimit

**Purpose:** Control the operating range of the tower crane system.

**Synopsis:** *Limit = get( tcrane, 'RailLimit' );*  
*set( tr3, 'RailLimit', NewLimit );*

**Description:** The *Limit* and *NewLimit* variables are 1x3 vectors. The elements of these vectors define the operating range of the cart and the maximum length of the lift-line. If a flag defined by the *RailLimitFlag* property is set the corresponding to it PWM wave stops when the corresponding to it encoder register exceeds the limit.

**See:** *RailLimitFlag*

## 8.8. RailLimitFlag

**Purpose:** Set range of limit flags.

**Synopsis:** *LimitFlag = get( tcrane, 'RailLimitFlag' );*  
*set( tr3, 'RailLimitFlag', NewLimitFlag );*

**Description:** The *RailLimitFlags* is a 1x3 vector. The first two elements control the operating range of the cart (a length of the rail and an angle position of the tower). The last element controls the maximum length of the lift-line. If the flag is set to 1 and the encoder register exceeds the range the DC motor corresponding to it stops. If the flag is set to 0 the motion continues in spite of the range limit exceeded in the encoder register.

**See:** *RailLimit, RailLimitSwitch*



## 8.9. RailLimitSwitch

**Purpose:** Read the state of limit switches.

**Synopsis:** *LimitSwitch = get( tcrane, 'RailLimitSwitch' );*

**Description:** Reads the state of three limit switches. Returns a 1x3 vector. If an element of this vector is equal to 0 it means that the switch has been pressed.

**See:** *RailLimit, RailLimitFlag*

## 8.10. ResetSwitchFlag

**Purpose:** Control the auto-reset of the encoder registers.

**Synopsis:** *ResetSwitchFlag = get( tcrane, 'ResetSwitchFlag' );*  
*set( tcrane, 'ResetSwitchFlag', ResetSwitchFlag );*

**Description:** The *ResetSwitchFlag* and *NewResetSwitchFlag* are 1x3 vectors. If an element of these vectors is equal to 1 the corresponding to it encoder register is automatically reset in the case when the corresponding to it limit switch is pressed.

**See:** *ResetEncoder, RailLimitSwitch*

## 8.11. Therm

**Purpose:** Read thermal flags of the power amplifiers.

**Synopsis:** *Therm = get( tcrane, 'Therm' );*

**Description:** Returns three thermal flags of three power amplifiers. When the temperature of a power amplifier is too high the appropriate flag is set to x.

**See:** *ThermFlag*

## 8.12. ThermFlag

**Purpose:** Control an automatic power down of the power amplifiers.

**Synopsis:** *ThermFlag = get( tcrane, 'ThermFlag' );*  
*set( tcrane, 'ThermFlag', NewThermFlag );*

**Description:** The *ThermFlag* and *NewThermFlag* are 1x3 vectors. If an element of these vectors is equal to 1 the DC motor corresponding to it is not excited by the PWM wave when it is overheated.

**See:** *Therm*

## 8.13. Time

**Purpose:** Return time information.

**Synopsis:**  $T = \text{get}(tcrane, 'Time');$

**Description:** The *towercrane* object contains the time counter. When a *towercrane* object is created the time counter is set to zero. Each reference to the *Time* property updates its value. The value is equal to the number of milliseconds past since the object was created.

#### 8.14. Quick reference table

Property Name	Description
<i>BitstreamVersion</i>	Read the version of the logic design for the RT-DAC/USB board
<i>Encoder</i>	Read the incremental encoder registers
<i>PWM</i>	Set the parameters of the PWM waves
<i>PWMPrescaler</i>	Determine the frequency of the PWM waves
<i>ResetEncoder</i>	Reset the encoder counters
<i>RailLimit</i>	Control the operating range of the Tower Crane system
<i>RailLimitFlag</i>	Set the range limit flags
<i>RailLimitSwitch</i>	Read the state of the limit switches
<i>ResetSwitchFlag</i>	Control the auto-reset of the encoder registers
<i>Therm</i>	Read the thermal flags of the power amplifiers
<i>ThermFlag</i>	Control the automatic power down of the power amplifiers
<i>Time</i>	Return time information

## 9. How to fill in the compilation settings page

The *Simulation/Configuration Parameters* page decides how works the process of compilation and linking. This page depends on the used third party compiler.

In Fig. 9.69 the *Configuration Parameters* page for the Visual C++ compiler is shown.




Fig. 9.69 *Configuration Parameters* page for the Visual C++ compiler

Notice that system target file *towercrane.tlc* is chosen. Also note that template makefile, responsible for compilation and linking process, is *towercrane\_vc.tmf*. It means the MS VC++ compiler is used in this case.

If the OpenWatcom compiler is used the *Configuration Parameters* page has different parameters. These parameters are shown in Fig. 9.70. Notice that system target file *towercrane.tlc* is the same as in the previous example. But the template makefile, responsible for compilation and linking process, is *towercrane\_openwatcom.tmf*

Interface page, responsible for connection with external target system, is shown in Fig. 9.71. To fulfil it properly a user must select *External mode* in *Interface* edit window. After that the user have to set *RT-CON tcpip* in the *Transport layer* edit window.




Fig. 9.70 Configuration Parameters page for the OpenWatcom compiler




Fig. 9.71 Interface page

In the *RT-CON options* a user can mark the *Display start-up message* checkbox (Fig. 9.72).




Fig. 9.72 RT-CON options window

If this checkbox is marked the message on loading real-time code is displayed (Fig. 9.73) after connecting to target in the model window.




Fig. 9.73 Start-up message